


## ramé-hart Model 190

**Contact Angle (CA) Goniometer** with **DROPimage CA Software** 

## ramé-hart Model 190 (p/n 190-U4)

Model 190 is an affordable and easy-to-use contact angle tool designed to offer the highest level of precision computer-driven contact angle analysis at the lowest possible acquisition cost. Model 190 accepts the Manual Tilting Base, Environmental Fixture, Wafer Support up to 4" in diameter, and many other options. This turn-key system includes DROPimage CA software.

What's in the box: Goniometer with U4 Series Camera (520 fps), Microsyringe Fixture and Shade, Halogen Fiber Optic Illuminator, (1) Microsyringe Assembly and (1) 22g Straight Needle for manual dispensing, PC with SuperSpeed Interface and monitor, DROPimage CA Software, User Manual, Calibration Tool, and Storage Cover — everything needed to start taking contact angle measurements right out of the box. Model 190 is a complete turn-key contact angle solution.



## **Specifications**

Stage Size		2 x 3 in (51 x 76mm)	
Sample Size		3 in (76mm) deep x unlimited	
Contact Angle Range		0 to 180°	
Resolution		0.1°	
Accuracy		<b>+/-</b> 0.10°	
U4 Series Camera		SuperSpeed USB 3.0, 1/3" CCD,	
		520fps, Sony S	ensor, 720 x 540 px
Backlighting		Variable Fiber Optic Illuminator	
Stage		Precision 3-Axis	s Locking Leveling
Dimensions		19 x 20 x 10 in	(480x500x250mm)
Weight	14 lbs / 6.3 kg (		excludes PC)
Power Supply		110 or 220 VAC	
DROPimage CA Tools Co		Contact Angle Measurement	
		Calibration Too	I I
Spares	Micros	yringe Assy.	100-10-20
	Needle		100-10-12-22
	Set of S	Spares	100-SPARES-F
	Fiber Optic Bulb		100-00-FOB
Options	Automated Dispensing		100-22
	Environmental Fixture		100-07-60
	Software Upgrade		400-22-30
	Manual Tilting Base		100-25-M
	Heated Enviro Cell		100-33
	Hot Plate		100-33-HP
	Film Clamps		100-15
	NIST Calibration Device 100-27-31-UPG		

ramé-hart instrument co• www.ramehart.com• carl@ramehart.com• 973-448-0305• fax 0315 19 Route 10 East•Suite 11•Succasunna•New Jersey•07876•USA